


JY/F-Q1 IU Balloon Tamponade Model

The model with accurate anatomical structure, mounted on a base, shows uterus outline on it. Comes with a bulb which can be inflated by water injection to simulate the postpartum uterus changes from 10 minutes to 48hours after placenta delivery. The access of the uterus cavity can be seen.



3 options of colors are available: Caucasian, Asian and African.



Minimum order quantity for your choice of color is 20 units.

Packing List	
Item Specification	Quantity
Soft portable bag	1
Instruction manual(in English,French and Spanish)	1
Warranty card	1